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WORK DETAILS

Engineer Name: MASASHI HIKASA

Event Type: Predictive / Preventative Maintenance

Device Status: Pass

Service Activity: Visit 2

EQUIPMENT DETAILS

Product Name: Ingenia 3.0T

Serial Number: 42263

Product Number: 781377

Installed Product: 57986182

UDI:

	Min	Max	UoM	Measurement	Result	Comments
PM section 4. Safety						
4. Safety						
Prevent electro static discharge					PASS	
PM section 6. Workflow						
On-site preparations						
Interview the representative of the site					YES	
PM section 7. System level						
System level						
Do PIQT (Remote / On-site)					PASS	
Check system parts visually					PASS	
Clean the computers externally					PASS	
Clean the EWS (Extended Workspace)					PASS	
Check the dust filter					PASS	
Clean the dust filter					PASS	
Check the dust filter of the patient ventilation					PASS	
PM - Disable alerts					PASS	
Check the RF enclosure					PASS	
Check the shock indicator(s)					PASS	
Check the RF coils for damage					PASS	
Relative humidity Examination room	40	70	%	44	Pass	

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

PM section 7. System level

System level

Relative humidity Technical room	30	70	%	44	Pass	
Temperature Examination room (> R5)	18	22	°C	21	Pass	
Temperature Technical room	15	24	°C	22	Pass	
Check the respiratory sensor and hose					PASS	
Check the version of the SPT spec file (Remote / On-site)					PASS	
Compare and download the SPT spec file					PASS	

PM section 9. Gradient System

Gradient System

Check the circuit to detect overheated gradient coils					PASS	
---	--	--	--	--	------	--

PM section 10. Liquid Cooling

Liquid cooling

Check the flow of the wc RF amplifier (WDH1.5T, AN8137, QLH_128)	16	20	l/minute	19	Pass	
Check the filling pressure	1.4	2.4	bar	2	Pass	
Check the hoses with liquid cooling for leaks					PASS	

Liquid cooling [LCC2BB] [LCC2B]

Fill the gradient coil and the WC RF amp (if < 1.4 bar)	1.4	2.4	bar	2	Pass	
Fill the gradient amp (if < 1.4 bar)	1.4	2.4	bar	2	Pass	
Check the flow in the gradient coil circuit	18	22	l/minute	20	Pass	
Check the flow in the gradient amplifier (C78X)	30	45	l/minute	40	Pass	

PM section 11. Refrigerator system

Refrigerator system

Check the compressor pressure F50 (dynamic)	1.9	2.2	MPa	2.1	Pass	
Check the water flow					PASS	
Check the cryogenic performance					PASS	

PM section 12. Magnet system

Magnet system

Check the labels and covers of the ERDU buttons					PASS	
Check the outlet of the helium venting system					PASS	
Check the service history of the coldhead (48 months)					≤ 48 Months Ago	
Record the last date coldhead service was done					Jul 29, 2020	

	Min	Max	UoM	Measurement	Result	Comments
PM section 12. Magnet system						
Magnet system						
Check Differential Pressure Switch (DPS)					PASS	
PM section 13. Patient support						
Patient support [IMT]						
Check the patient alarm					PASS	
Check the patient comfort lights					PASS	
Clean the DCI RF coil sockets and connector					PASS	
Clean the PSU					PASS	
PM section 14. System Meas. & Adjustments						
System Meas. & Adjustments						
Check the body coil transmit characteristics					PASS	
Run the software tests for PM						
System selftest					PASS	
PFEI peripheral test					PASS	
Patup self test					PASS	
PSU2 selftest (PSU(2) fan speed test)					PASS	
TX att. calibration					PASS	
Max kW calibration body coil					PASS	
RF amp TX1 PMU test body coil					PASS	
RF amp TX2 PMU test body coil					PASS	
RF amp controlloop calibration					PASS	
F0 determination					PASS	
RF power ref. cal. body coil					PASS	
Pickup coil tripl.cal.body coil					PASS	
Multitransmit RF calibration					PASS	
System measurements and adjustments						
Check the VSWR (3.0T multi-transmit TX1)		5		1.4	Pass	
Check the VSWR (3.0T multi-transmit TX2)		5		1.4	Pass	
Use the analyzer					PASS	
Run the RF spurious test					PASS	
Run the spurious noise test batch					PASS	
Analyze the result of the spurious noise test batch					PASS	
PIQT						
Do PIQT					PASS	

	Min	Max	UoM	Measurement	Result	Comments
PM section 14. System Meas. & Adjustments						
PIQT						
Generate a report with PIQT results					PASS	
PM section 15. Operator Console & Viewing						
Operator Console & Viewing (All)						
Clean the mouse, keyboard, monitor of host and EWS					PASS	
Check the audio communications					PASS	
PM section 16. Software						
Software						
Check and update the EAS software (e-Alert)					PASS	
Clean up disks					PASS	
Defragment the hard disks					PASS	
Do the PRS related checks					PASS	
Make a system back-up					PASS	
Export the configuration					PASS	
PM section 19. PM end						
19. PM end						
PM - Enable alerts					PASS	
Handover					PASS	
Administration					PASS	

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