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Chapter: Liquid Cooling Eaton Williams (All)

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Chapter: Liquid Cooling Eaton Williams (1, 3)

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Chapter: Liquid Cooling Eaton Williams (1)

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Chapter: Operator Console & Viewing (All)

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Chapter: Workflow (All)

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Chapter: System Level Checks (Visits 1,3)

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WORK DETAILS

Engineer Name: YUSUKE SADAOKA

Event Type: Predictive / Preventative Maintenance

Device Status: Pass

Service Activity: Additional Visit

EQUIPMENT DETAILS

Product Name: Achieva 1.5T

Serial Number: 32625

Product Number: 781296

Installed Product: 51640158

UDI:

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

PM section 4. Safety

4. Safety

Prevent electro static discharge					PASS	
----------------------------------	--	--	--	--	------	--

PM section 6. Workflow

Workflow

Interview the representative of the site					YES	
--	--	--	--	--	-----	--

PM section 7. System level

System level

Do PIQT (Remote / On-site)					PASS	
Check system parts visually					PASS	
Check the dust filter					PASS	
Check the dust filter of the patient ventilation					PASS	
Disable alerts					PASS	
Check the RF enclosure					PASS	
Check the RF coils for damage					PASS	
Relative humidity Examination room	40	70	%	65.4	Pass	
Relative humidity Technical room	30	70	%	54.5	Pass	
Temperature Examination room (R3 or upgraded from R3)	20	24	°C	22.3	Pass	
Temperature Technical room	15	24	°C	21.4	Pass	
Check the version of the SPT spec file (Remote / On-site)					PASS	

	Min	Max	UoM	Measurement	Result	Comments
PM section 7. System level						
System level						
Compare and download the SPT spec file					PASS	
PM section 10. Liquid Cooling						
Liquid cooling						
Check the flow in the gradient coil circuit	18	22	l/min	22	Pass	
Check the flow in the gradient amplifier	30	45	l/min	34	Pass	
Liquid cooling [LCC2]						
Check the hoses with liquid cooling for leaks					PASS	
PM section 11. Refrigerator system						
Refrigerator system						
Check the dynamic pressure of the compressor					PASS	
PM section 12. Magnet system						
Magnet system						
Check the labels and covers of the ERDU buttons					PASS	
Check the service history of the coldhead [10K]					≤ 36 Months Ago	
Record the last date coldhead service was done					Jan 11, 2023	
PM section 13. Patient support						
Patient support [MT]						
Check the stop function on operator console					PASS	
Check the stop scan buttons on the magnet					PASS	
Check the emergency stop circuit					PASS	
Check the stop function of finger protection plate					PASS	
Check the patient alarm					PASS	
Check the patient comfort lights					PASS	
PM section 14. System Meas. & Adjustments						
System Meas. & Adjustments						
Do a check with the HealthViewer application					Not applicable (< R5.3.1 and/or CDAS)	
System measurements and adjustments						
Use the analyzer					PASS	

Min Max UoM Measurement Result Comments

PM section 14. System Meas. & Adjustments

PIQT

Do PIQT	PASS
Generate a report with PIQT results	PASS
Mark the PIQT scans as reference specifications	PASS

PM section 15. Operator Console & Viewing

Operator Console & Viewing (All)

Clean the mouse, keyboard, monitor of host and EWS	PASS
Check the audio communications	PASS

PM section 16. Software

Software

Check and update the EAS software (e-Alert)	PASS
Clean up disks	PASS
Defragment the hard disks	PASS
Export the configuration	PASS

PM section 17. FCOs and corrective maintenance

FCOs and corrective maintenance

Implement FCOs and corrective maintenance	N/A
---	-----

Philips Representative Signature:

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