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WORK DETAILS

Engineer Name: MASASHI HIKASA

Event Type: Predictive / Preventative Maintenance

Device Status: Pass

Service Activity: Visit 3

EQUIPMENT DETAILS

Product Name: Achieva 1.5T

Serial Number: 32625

Product Number: 781296

Installed Product: 51640158

UDI:

	Min	Max	UoM	Measurement	Result	Comments
PM section 4. Safety						
4. Safety						
Prevent electro static discharge					PASS	
PM section 6. Workflow						
Workflow						
Interview the representative of the site					YES	
PM section 7. System level						
System level						
Do PIQT (Remote / On-site)					PASS	
Check the door switch					Pass	
Check the system airflow [SACU]					Not Installed	
Check system parts visually					PASS	
Check the safety signs					PASS	
Check the dust filter					PASS	
Check the dust filter of the patient ventilation					PASS	
Disable alerts					PASS	
Check the RF enclosure					PASS	
Check the RF coils for damage					PASS	
Relative humidity Examination room	40	70	%	54	Pass	
Relative humidity Technical room	30	70	%	55	Pass	

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

PM section 7. System level

System level

Temperature Examination room (> R5)	18	22	°C	21	Pass	
Temperature Technical room	15	24	°C	21	Pass	
Check the respiratory sensor and hose					PASS	
Check the version of the SPT spec file (Remote / On-site)					PASS	
Compare and download the SPT spec file					PASS	

PM section 9. Gradient System

Gradient System

Check the connections of the gradient coil					PASS	
--------------------------------------------	--	--	--	--	------	--

PM section 10. Liquid Cooling

Liquid cooling

Check the flow in the gradient coil circuit	18	22	l/min	20	Pass	
Check the flow in the gradient amplifier	30	45	l/min	33	Pass	

Liquid cooling [LCC2]

Check the hoses with liquid cooling for leaks					PASS	
Flush the circuit of the gradient coil	2.2	2.4	bar	2.3	Pass	
Flush the gradient amp	2.2	2.4	bar	2.2	Pass	
Check the flow of the wc RF amplifier (WDH1.5T, AN8137)	16	20		20	Pass	
Check the flow of the wc RF amplifier (S35)	9	11		10	Pass	
Check the flow Interlock					PASS	

PM section 11. Refrigerator system

Refrigerator system

Check the dynamic pressure of the compressor					PASS	
Check the cryogenic performance					PASS	
Check for leaks					PASS	

PM section 12. Magnet system

Magnet system

Check the labels and covers of the ERDU buttons					PASS	
Check the water drain					PASS	
Check the bolts of the helium vent pipe					PASS	
Check the service history of the coldhead [10K]					< 24 Months	
Record the last date coldhead service was done					Oct 9, 2025	

Min Max UoM Measurement Result Comments

PM section 12. Magnet system

Magnet system

Record the compressor type				HC-8E	
Record the Coldhead serial number				A8243	

PM section 13. Patient support

Patient support [MT]

Check the stop function on operator console				PASS	
Check the stop scan buttons on the magnet				PASS	
Check the emergency stop circuit				PASS	
Check the emergency stop button				PASS	
Check the tabletop release button				PASS	
Check the stop function of finger protection plate				PASS	
Check the manual function of the tabletop				PASS	
Check the patient alarm				PASS	
Clean the pollution of grease				PASS	
Check the patient comfort lights				PASS	
Check the table top for damage, loose glass-fibers				PASS	
Check the wheels				PASS	
Check the patient support				PASS	

PM section 14. System Meas. & Adjustments

System measurements and adjustments

Use the analyzer				PASS	
------------------	--	--	--	------	--

PIQT

Do PIQT				PASS	
Generate a report with PIQT results				PASS	
Mark the PIQT scans as reference specifications				PASS	

PM section 15. Operator Console & Viewing

Operator Console & Viewing (All)

Clean the mouse, keyboard, monitor of host (and ISP/EWS)				PASS	
Check the audio communications				PASS	
Maintain the ceiling carriage				PASS	

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

PM section 16. Software

Software

Check and update the EAS software (e-Alert): check the version of the software on the system				e-Alert is not installed	
Find the latest version of the EAS software (e-Alert) on InCenter				https://philips.mizecx.com/knowledgeDocument.html?17eeaa963b1881a0c0354f5a426d4f93	
Clean up disks				PASS	
Defragment the hard disks				PASS	
Export the configuration				PASS	

PM section 17. FCOs and corrective maintenance

FCOs and corrective maintenance

Implement FCOs and corrective maintenance				FCO78100610	
-------------------------------------------	--	--	--	-------------	--

PM section 19. PM end

19. PM end

PM - Enable alerts				PASS	
Handover				PASS	
Administration				PASS	

Philips Representative Signature:

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