

RESULTS OF OSLD CHECK OF PHOTON BEAM OUTPUT

v 8.0.2


Institution:	Shikoku Cancer Center, Matsuyama City, Ehime Prefecture
RTF Number	4366
Person irradiating dosimeters:	Kazunobu Koshi
Radiation Machine:	Clinac 21EX Serial 3132
Radiation Quality:	4 MV X-rays
Distance from source to reference point:	100.0 cm

OUTPUT VERIFICATION:

Date of Irradiation	IROC Houston measured dose at dmax:*	Institution reported dose at dmax:*	Ratio of absorbed dose determined by IROC Houston to that stated by institution: OSLD/INST
12-Oct-2017	99.4 cGy to water	99.9 cGy to water	1.00

Agreement within 5% is considered a satisfactory check.

RESULT HISTORY FOR THIS BEAM



THIS INFORMATION SHOULD BE USED ONLY AS A CHECK OF MACHINE OPERATION AND NOT AS A MACHINE CALIBRATION, nor as an alternative to frequent calibration by a qualified physicist.


The OSLD dose was evaluated using the AAPM TG-51 Dosimetry Calibration Protocol.

OSLD read on: 27-Oct-2017
 OSLD read by: Omar Garcia
 Checked by: Stephen Kry, Ph.D.

*The variance of the dose determined by a single OSLD is less than 3%. The OSLD sample, therefore, has an uncertainty of 5% at a confidence level in excess of 90%. This analysis did not include uncertainties in the institutions' irradiation technique.

David S. Followill

 David S. Followill
 Director



IROC Houston QA Center
 MD Anderson Cancer Center
 8060 El Rio Street
 Houston, TX 77054
 Tel (713) 745-8989
 Fax (713) 794-1364
 Email: irochouston@mdanderson.org

RESULTS OF OSLD CHECK OF PHOTON BEAM OUTPUT
 v 8.0.2


Institution:	Shikoku Cancer Center, Matsuyama City, Ehime Prefecture
RTF Number	4366
Person irradiating dosimeters:	Kazunobu Koshi
Radiation Machine:	Clinac 21EX Serial 3133
Radiation Quality:	4 MV X-rays
Distance from source to reference point:	100.0 cm

OUTPUT VERIFICATION:

Date of Irradiation	IROC Houston measured dose at dmax:*	Institution reported dose at dmax:*	Ratio of absorbed dose determined by IROC Houston to that stated by institution: OSLD/INST
12-Oct-2017	99.3 cGy to water	99.9 cGy to water	0.99

Agreement within 5% is considered a satisfactory check.

RESULT HISTORY FOR THIS BEAM



THIS INFORMATION SHOULD BE USED ONLY AS A CHECK OF MACHINE OPERATION AND NOT AS A MACHINE CALIBRATION, nor as an alternative to frequent calibration by a qualified physicist.

The OSLD dose was evaluated using the AAPM TG-51 Dosimetry Calibration Protocol.


OSLD read on: 27-Oct-2017
 OSLD read by: Omar Garcia
 Checked by: Stephen Kry, Ph.D.


David S. Followill
 Director

*The variance of the dose determined by a single OSLD is less than 3%. The OSLD sample, therefore, has an uncertainty of 5% at a confidence level in excess of 90%. This analysis did not include uncertainties in the institutions' irradiation technique.



IROC Quality Assurance Center Locations
 Houston | Ohio | Philadelphia | Rhode Island | St. Louis
 Sponsored by the National Cancer Institute





IROC Houston QA Center
 MD Anderson Cancer Center
 8060 El Rio Street
 Houston, TX 77054
 Tel (713) 745-8989
 Fax (713) 794-1364
 Email: irochouston@mdanderson.org

RESULTS OF OSLD CHECK OF PHOTON BEAM OUTPUT
 v 8.0.2


Institution:	Shikoku Cancer Center, Matsuyama City, Ehime Prefecture
RTF Number	4366
Person irradiating dosimeters:	Kazunobu Koshi
Radiation Machine:	Clinac 21EX Serial 3133
Radiation Quality:	10 MV X-rays
Distance from source to reference point:	100.0 cm

OUTPUT VERIFICATION:

Date of Irradiation	IROC Houston measured dose at dmax:*	Institution reported dose at dmax:*	Ratio of absorbed dose determined by IROC Houston to that stated by institution: OSLD/INST
12-Oct-2017	100.8 cGy to water	100.0 cGy to water	1.01

Agreement within 5% is considered a satisfactory check.

RESULT HISTORY FOR THIS BEAM



THIS INFORMATION SHOULD BE USED ONLY AS A CHECK OF MACHINE OPERATION AND NOT AS A MACHINE CALIBRATION, nor as an alternative to frequent calibration by a qualified physicist.

The OSLD dose was evaluated using the AAPM TG-51 Dosimetry Calibration Protocol.

OSLD read on: 27-Oct-2017
 OSLD read by: Omar Garcia
 Checked by: Stephen Kry, Ph.D.


David S. Followill
 Director

*The variance of the dose determined by a single OSLD is less than 3%. The OSLD sample, therefore, has an uncertainty of 5% at a confidence level in excess of 90%. This analysis did not include uncertainties in the institutions' irradiation technique.



IROC Quality Assurance Center Locations
 Houston | Ohio | Philadelphia | Rhode Island | St. Louis
 Sponsored by the National Cancer Institute





IROC Houston QA Center
 MD Anderson Cancer Center
 8060 El Rio Street
 Houston, TX 77054
 Tel (713) 745-8989
 Fax (713) 794-1364
 Email: irochouston@mdanderson.org

RESULTS OF OSLD CHECK OF PHOTON BEAM OUTPUT

v 8.0.2


Institution: Shikoku Cancer Center, Matsuyama City, Ehime Prefecture
RTF Number: 4366
Person irradiating dosimeters: Kazunobu Koshi
Radiation Machine: Clinac 21EX Serial 3132
Radiation Quality: 10 MV X-rays
Distance from source to reference point: 100.0 cm

OUTPUT VERIFICATION:

Date of Irradiation	IROC Houston measured dose at dmax:*	Institution reported dose at dmax:*	Ratio of absorbed dose determined by IROC Houston to that stated by institution: OSLD/INST
12-Oct-2017	100.7 cGy to water	100.0 cGy to water	1.01

Agreement within 5% is considered a satisfactory check.

RESULT HISTORY FOR THIS BEAM



THIS INFORMATION SHOULD BE USED ONLY AS A CHECK OF MACHINE OPERATION AND NOT AS A MACHINE CALIBRATION, nor as an alternative to frequent calibration by a qualified physicist.

The OSLD dose was evaluated using the AAPM TG-51 Dosimetry Calibration Protocol.

OSLD read on: 27-Oct-2017
 OSLD read by: Omar Garcia
 Checked by: Stephen Kry, Ph.D.


David S. Followill
 Director

*The variance of the dose determined by a single OSLD is less than 3%. The OSLD sample, therefore, has an uncertainty of 5% at a confidence level in excess of 90%. This analysis did not include uncertainties in the institutions' irradiation technique.



IROC Quality Assurance Center Locations
 Houston | Ohio | Philadelphia | Rhode Island | St. Louis
 Sponsored by the National Cancer Institute





The information contained in this report may be privileged, confidential, and protected from disclosure. If you are not the intended recipient, any further disclosure or use, dissemination, distribution, or copying of this report is strictly prohibited. If you think that you have received this report in error, please delete it and notify the sender.