

2025年 6 月

検査機器保守管理作業日誌（肺機能1）

部門責任者

		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
始業時	外観の異常がないこと		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	コード類の破損がないこと		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	電源投入時の表示が通常通り		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	機器の動作確認		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	肺機能環境データの確認		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	メッシュセンサの装着		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	肺機能キャリブレーション （ファクター：0.90～1.10以内）		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	肺機能精度管理 （測定データ：±3%以内）		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	バイオミクロンサークル（空気清浄機）の 外観に異常がない・電源ON		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	時間		8:30	8:50	8:20	8:30	8:30			8:30	8:30	8:20	8:35	8:20			8:20	8:20	8:30	8:20	8:20			8:30	8:20	8:25	8:10	8:30			8:30	
確認印	[Redacted]																															

終了時	使用済みスパイロフィルタ等の破棄，密閉		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	メッシュセンサの破棄		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	バイオミクロンサークル（空気清浄機） 電源OFF		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	消耗品の在庫確認		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	
	時間		17:30	17:50	17:35	18:15	18:30			17:30	17:35	17:25	17:45	17:50			17:30	17:30	18:15	17:30	17:30			17:35	17:40	18:15	17:45	17:30			17:40	
確認印	[Redacted]																															

備考	修理・メンテナンス 機器清拭/動作確認	[Redacted]																														
	バイオミクロンサークル（空気清浄機）の フィルター等の清掃（月1回）	[Redacted]																														
		[Redacted]																														

✓：問題なし ×：確認を行う