

部門責任者

2025年6月

検査機器保守管理作業日誌（超音波6）

EPIQ Elite		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31		
始業時	外観の異常がないこと		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓			
	コード、プローブ類の破損がないこと		✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓			
	電源投入時の表示が通常通り		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	機器の動作確認		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	電動ベッドの動作確認		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	心電図ディスプレイ装着		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	ゲルウォーマー電源ON		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	時間		8:20	8:20	8:25	8:20	8:15				8:20	8:20	8:20	8:20	8:20			8:25	8:20	8:20	8:25	8:15			8:20	8:15	8:15	8:20	8:15			8:20		
確認印																																		
終了時	消耗品の在庫確認		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	血液、体液による汚染がない		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	ゲルウォーマーの電源OFF		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	プローブの清拭		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	心電図ディスプレイ取り外し		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
	おしぼりの回収		✓	✓	✓	✓	✓				✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓	✓	✓	✓	✓			✓		
時間		17:30	17:30	17:25	18:10	17:20				17:35	17:25	17:15	18:00	18:10			17:25	17:35	17:50	17:30	17:15			17:35	17:35	17:15	17:30	17:35			17:40			
確認印																																		
その他	枕カバー、タオルの交換		16:25	17:40	17:10	17:40	17:20				17:15	16:55	17:15	17:20	17:05			15:25	16:55	16:55	16:50	16:00			17:10	17:35	15:55	17:30	16:10			17:15		
	超音波装置、フィルターの清掃																																	
	心電図電極、プローブホルダーの洗浄																																	
	修理・メンテナンス 機器清拭/動作確認																																	
備考																																		
※時刻合わせを行った翌日の始業点検時には時刻を確認		時刻合わせOK																																